# Lovibond® Water Testing

Tintometer® Group

## Floc-Tester

### ET 730

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Seite 3–12

**GB** Instruction Manual

Page 13–22

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Page 23–32

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Pagina 33–42

[www.lovibond.com](http://www.lovibond.com)

---

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---

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- **General Notes** ..... 14
  - Safety Instructions..... 15
  - Introduction ..... 15
  - Instrument description ..... 15
  
- **Functional description** ..... 16
  
- **Assembly and installation** ..... 17
  - Electric connection ..... 17
  - Start-up ..... 17
  
- **Operating controls** ..... 18
  
- **User notes** ..... 19
  
- **Calculation of velocity gradient (G)** ..... 20
  
- **Technical data** ..... 21
  
- **Declaration of conformity** ..... 22

## Safety Instructions

1. Before using the unit, please read the operating manual supplied with the apparatus carefully.
2. Do not dispose of this equipment as normal household waste.
3. In order to prevent possible risk of electric shocks, fire and personal injury when the unit is being used, basic safety measures must always be taken, including:
4. Do not turn the instrument on before positioning the beakers (thus protecting the blades). Ensure that the power switch is turned to "0" (OFF) before removing the beakers containing the stirring liquid.
5. Check that the rating of external power supply corresponds to the rating of the electric line.
6. Do not use the unit if it is not working correctly. In case of malfunctioning, contact the manufacturer.
7. Personal protective equipment must be compatible with the possible risks posed by the material being processed and the glass parts.
8. Follow the cleaning instructions described in this manual.
9. This unit must only be used for laboratory applications.
10. The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

## Note

- The manufacturer is committed to constantly improving the quality of the products and reserves the right to modify the characteristics without prior notice.

## Introduction

Jar Testers are studied both to optimize the dosing of coagulants in order to separate the pollutions in water treatment plants, thanks to the laboratory valuation results obtained with the so-called Jar Test.

It is possible to set the stirring speed in analogical way and the stirring time by knobs placed on the frontal panel.

The examined samples can be lighted by a back fluorescent lamp controlled by a switch located on front panel.

The portable jar tester ET 730 was developed for optimum use in the field, in fact the instrument can run both connected to a normal battery and a car's cigarette lighter.

## Instrument description

The steel structure epoxy painted is studied and tested in order to give to the instrument an excellent resistance to the attack of chemical agents and to corrosion in general.

The external power supply allows the instrument to be connected to all frequencies and voltages in the world (100÷240 V – 50÷60 Hz).

It's possible to select the stirring speed (20-40-50-100-200 rpm), the time of running and the lighting of central column by the frontal panel. The stirring speed is the same for each position. It is possible to set the time of running from 0 to 30 minutes and the central column start up.

ET 730 is characterized by a non-skid base so as to keep the beakers stable during operation.

---

## **GB** Functional description

---

### **Stirring Speed Knob**

Turn the stirring speed knob on the front panel in order to start the stirring function. It is possible to select the following speeds: 0 – 20 – 40 – 50 – 100 – 200 rpm

### **Timer**

It's possible to select the time of running by the Timer knob on the frontal panel: turn counter clockwise for continuous running or turn clockwise in order to select the time up to 30 minutes. At the end of time the instrument stops automatically the running.

### **Lighting Column Switch**

It's possible to turn on the light of central column by the switch.

### **Information regarding construction materials**

Structure	Steel
Frontal panel label	PET
Foot	NBR

If any liquid being processed should fall onto the above materials, clean them immediately to avoid chemical corrosion.

## **Electric connection**

After unpacking the instrument, place the unit on a laboratory bench.

Before connecting the instrument to the electric line, make sure that the values on the rating plate correspond to those of the electric line.

Use the external power supply to connect the instrument to the socket.

The power supply can be from 100 to 240 V AC with a frequency of 50 or 60 Hz.

## **Start-up**

After the connection of the instrument to the electric feeding place the beakers with sample on the instrument with inside the stirring shaft.

Pull to the top the stirring shaft by the black plastic knob in order to insert the beakers in each position and then take them down.

Do not turn the instrument on before positioning the beakers (thus protecting the blades). Ensure that the power switch is turned to "0" (OFF) before removing the beakers containing the stirring liquid.

---

## **GB** Operating controls

---

It's possible to select the desired stirrer speed by the knob "Stirrer rpm" on the frontal panel. It's possible to chose the desired values: 0-20-40-50-100-200 rpm by using the information written on the analogical scale around the knob.

The start up and the selection of the running time with automatical turn off are made by the relative timer knob:

by turning unclockwise the timer knob on the manual modality the instrument runs in continuous way;

by turning clockwise the timer knob it's possible to select the time of running up to 30 minutes; the stirring automatically stops as the time finishes.

## **End of work cycle operations**

At the end of the work cycle set the stirring speed knob on "0" rpm, the timer knob on "OFF" and turn off the light of central column.

If the instrument is not used for a long period it's better to disconnect the external power supply from the electric line.

## **Maintenance**

No routine or extraordinary maintenance is necessary but periodically cleaning of the unit as described in this manual.

In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

## **Cleaning**

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable and non-aggressive detergent to clean the unit.

## **Disposal**

The final disposal of the unit or of its components must be carried out in compliance with the directives and laws in force in the country concerned.

## GB Calculation of velocity gradient (G)

### Formula

Velocity gradient ( $s^{-1}$ )

$$G = \sqrt{P/\mu V}$$

Power requirement (W)

$$P = k * \rho * n^3 * D^5 \text{ turbulent}$$

$$P = k * \mu * n^2 * D^3 \text{ laminar}$$

### Notes


Dynamic viscosity:	$\mu$	0,001002	( $N*s*m^{-2}$ )
Mass density of fluid:	$\rho$	998,2	( $kg*m^{-3}$ )
Water volume:	V	0,001	( $m^3$ )
Liquid mixed	Water		
Measurement temperature:	20		( $^{\circ}C$ )
Constant	k	3,38	
Number of turns / minute (RPM)	n	from 30 to 200	( $min^{-1}$ )
Number of turns / second	n	(from 30 to 200)*1/60	( $s^{-1}$ )
Diameter of impeller	D	0,075	(m)

### Example: Velocity gradient for turbulent flow

#### Data

Number of turns / min-1	G (turbulent) / $s^{-1}$
200	544,28
150	353,52
100	192,43
75	124,99
50	68,03
40	48,68
30	31,62

#### Graph

---

**GB** Technical data

---

**General**

Power supply	V	DC 12 ± 0,5
Power	W	6
Dimensions	mm (bxhxp)	250 x 305 x 250
Weigh	Kg	3,8
Programmable speed	rpm (1/min)	0 – 20 – 40 – 50 – 100 - 200
Construction material		Stell
Environmental temp. range	°C	+ 5...+ 40
Temperatura di stoccaggio ammessa	°C	- 10...+ 60
Max. humidity	%	Max 80
Operating mode		Continuous
Pollution degree CEI EN61010-1		2
Installation category CEI EN61010-1		2

**External Power Supply**

Inlet	AC 100±240 V ; 50-60 Hz ; 1.5 A
Outlet	DC 12 V ; 1.25 A

---

**GB Declaration of conformity**

---

The manufacturer    Tintometer GmbH  
Adress                Schleefstraße 8-12  
                             44287 Dortmund  
                             Germany

under our responsibility declare that the product is manufactured in conformity with the following standards:

EN 61010-1 (2001)

EN 61326-1 (1997) + A1 (1998) + A2 (2001) + A3 (2003)

2002/95/CE (RoHS)

2002/96/CE (WEEE)

and satisfies the essential requirements of the following directives:

Machines directive 2006/42/EC

Low voltage directive 2006/95/EC

Electromagnetic compatibility directive 2004/108/EC

plus modifications and that the documents listed in annex V are available as foreseen by the machine directive.

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**Tintometer Indien Pvt. Ltd.**