DE

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GB

Test Instructions

1. Use the photometer as described in the instructions for the method which correspond with normally used reagent system (tablet, liquid reagent or powder pack). For liquid reagent use the tablet mode.
2. Clean the outside of the vial with a clean tissue.
3. Zero the instrument according to instrument instructions using the "Reference Zero".
4. Use the Reference Standard instead of the test sample and take the test reading.
The reading should correspond with the stated concentration in relation to the used reagent system within the given tolerances (standard & instrument).

Notes

1. If an instrument and the Reference Standards show Non-Compliance:
 - the cell chamber should be checked for dirtiness, if so, cleaned and the measurement should be repeated
 - the instrument and the reference Standards should be returned to Tintometer GmbH, Technical Service Department, or an authorized distributor, for servicing
2. Place standards in the case after use.
Store in the dark at room temperature (15 - 25°C).
Do not use the standards after the stated expiry date.
3. Reference Standards should be used under laboratory conditions.
Avoid working in direct sunlight and take the readings at a temperature of $20 \pm 3^\circ\text{C}$.