

Lovibond® Water Testing Tintometer® Group



DI20 Incubator

Compliant with WHO guidelines



Our new DI20 incubator has been designed to allow the user to reach global drinking water standards for testing, even in difficult environments.

The DI20 incubator is a high quality, portable, lightweight incubator for microbiological testing of water using plates. It is the only incubator of its type to come with heating and cooling as standard ensuring the samples are incubated at the correct temperature no matter what ambient conditions.

Its innovative design also includes a small petri dish holder utilising our re-usable stainless steel micro plates, the dishes allow each incubator to perform up to 20 tests.

Its sleek design and intuitive pictographic buttons ensure that it can be used with little or no training and its LED lights allow users to see at a glance how far through a cycle the incubation is, and that the temperature is stable.

- Compliant with WHO guidelines " Guidelines for drinking water quality 4th edition."
- Compact design & lightweight
- Heating & Cooling fitted as standard
- Intuitive, unique pictographic buttons

Part Number: 56B000714

Industry

NGO

Technical Data

Interfaces	USB
Power Supply	12 V
Temperature Range	20 °C - 47 °C
Temperature Tolerance	± 0.5 °C
Environmental Conditions	5 °C - 50 °C
Compliance	CE
Dimensions	116 x 165 x 116 mm

Delivery Scope

- DI20
- Power leads
- 50 Petri dishes
- Power Supply

Title	Part Number
Power cable set	136300
USB Cable 2.0 A - Micro B plug	190630
Petri dish holder	197139
Power supply, 12 V / 48 W	19803550
Battery cable	19803555
Petri dishes, 10 pc.	400855

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germany

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
sales@lovibond.uk
www.lovibond.com
UK

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
USA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain