Lovibond® Colour Measurement

Tintometer® Group



Filterset IP-eenheden



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Filterset IP-eenheden

A set of four glass conformance filters for checking the PFXi instruments in IP (method B) units

Bestelnr.: 137660

Technische gegevens

 Glass Standard Values
 Approximate values

 es
 - 0.75

 - 1.25
 - 2.0

 - 3.5

Filters

Leveromvang

- ── IP 0.75 filter
- IP 1.25 filter
- ── IP 2.0 filter
- IP 3.5 filter
- Certificate of conformity