

Lovibond® Colour Measurement Tintometer® Group



Conjunto de filtros Unidades IP

Artigo No: 137660



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Conjunto de filtros Unidades IP

A set of four glass conformance filters for checking the PFXi instruments in IP (method B) units

Dados Técnicos

| | |
|------------------------------|------------------------|
| Glass Standard Values | Valores aproximados: - |
| | 0.75 |
| | 1.25 |
| | 2.0 |
| | 3.5 |

| | |
|----------------|---|
| Filters | 4 |
|----------------|---|

Texto do Escopo de Entrega

- IP 0.75 filter
- IP 1.25 filter
- IP 2.0 filter
- IP 3.5 filter
- Certificate of conformity