


# Lovibond® Colour Measurement

Tintometer® Group



## Daylight 2000 Unit

- ✓ Improved appearance and ergonomics
- ✓ Constructed from ABS (Acrylonitrile Butadiene Styrene) Fire Retardant Plastic. This robust, chemically resistant material is ideal for industrial environments
- ✓ Large switch illuminates sample only when pressed - no energy wastage
- ✓ Halogen, long life, light source combined with standard daylight filters ensures full compliance with former units and consistency from unit to unit
- ✓ Low voltage power supply

Daylight 2000 Lighting Unit : 171010  
Nessleriser Daylight 2000 : 171020