

Lovibond® Colour Measurement

Tintometer® Group



TRA 520 / TRA 500

- Portable docking station for use with the Lovibond® TR 520 / 500 series of spectrophotometers
- Provides a light controlled environment and sample optical cells for the colour measurement of liquids, gels, powders and other 'awkward' objects
- No ambient light - tight-fitting integral lid
- Supports a wide variety of samples with cell sizes of 10, 20 and 30mm
- Supplied with a white background plate to ensure consistency of readings
- Light weight and portable, yet very sturdy and robust
- Quick and easy to calibrate the instrument within the adapter
- Available as a package - TR 500 (TRA 500) and TR 520 (TRA 520) or as an upgrade from the TR 500 / 520 (serial number required)

Market sectors

Coatings

Food and Beverage

Household Products

Medical and Clinical

Plastics



TRA 520 / TRA 500 Innovative Product Design

The TRA 500 / 520 docking station is precision engineered to house the TR500 / TR 520 spectrophotometer at the optimum ergonomic position for touch screen operation and sample colour measurement.



A cover slots securely into place above the glass cell to ensure no ambient light can affect the reading



A white background plate is included to ensure consistency of colour measurement

10, 20 and 30mm optical glass cells are supplied with the unit to securely position and hold a variety of liquid, gel and powder samples

Order No	Product	Order No	Product
403225	Lovibond® TRA 520	605910	W586.OG.30mm Cell
403220	Lovibond® TRA 500 (8mm aperture)	605895	W586.OG.10mm Cell (Pack of 5)
605890	W586.OG.10mm Cell	605905	W586.OG.20mm Cell (Pack of 5)
605900	W586.OG.20mm Cell	605915	W586.OG.30mm Cell (Pack of 5)

www.lovibond.com

sales@tintometer.com

[Download form](#)

The Tintometer Ltd
T: +44 (0)1980 664800
F: +44 (0)1980 625412

Tintometer Inc
T: +1 941 756 6410
F: + 1941 727 9654

Tintometer GmbH
T: +49 (0) 231/94510-0
F: +49 (0) 231/94510-20

Tintometer – China
T: +86 10 8525 1111 Ext 330
F: +86 10 85251001

Tintometer – South East Asia
T: +603 3325 2286
F: +603 3325 2287

Tintometer India Pvt. Ltd
T: +91 (0) 40-23883300
Toll Free: 1800 599 3891/92

Lovibond® and Tintometer® are Registered Trade Marks of The Tintometer® Group. All translations and transliterations of LOVIBOND® & TINTOMETER® are asserted as Trade Marks of The Tintometer Limited. Registered Office: Lovibond House, UK. Registered in England No. 45024.

Errors and Omissions Excepted – Content Subject to alterations without notice.

The Bluetooth® mark, figure mark and combination mark are trademarks owned by the Bluetooth SIG and any use by the Tintometer® Group is under license. Please note, Bluetooth® is not available for all regions. Please contact the local Lovibond® representative for details.