

Lovibond® Water Testing

Tintometer® Group



PDM/MAC Dipslide, Coliforms/Pseudomonas, Dual Sided, 10 pack

With dual agar for the determination of Pseudomonas / Enterobacteriaceae



- Determine aerobic and anaerobic bacteria levels
- Early indication of bacteria proliferation
- Results in 48 hours
- Inexpensive

Part Number: 56B010610

Guidelines to advise on the correct practices to control legionella bacteria in water systems exist worldwide. Suppressing bacteria levels does substantially reduce the possibility of an outbreak of this often fatal disease. Guidelines recommend introducing a monitoring and control program. This program includes the testing of cooling tower waters with dipslides on a regular basis. A full range of dipslides is available for semiquantitative determination of aerobic and anaerobic bacteria populations in industrial and recreational waters. Dipslide accuracy is limited due to the small sample size, but if used correctly and incubated at a constant temperature using the Lovibond® dipslide incubator, they are excellent for trend analysis and can give an early indication of bacteria proliferation.

Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Marine Industry | Oil and Gas | Pharmaceutical Industry | Power and Energy

PDM/MAC Dipslide, Coliforms/Pseudomonas, Dual Sided, 10 pack

If you want to detect Pseudomonas and Enterobacteriaceae simultaneously, this dipslide is the right choice. The paddles are coated with an agar containing Pseudomonas base medium and MacConkey No. 3 agar. The dipslide corresponds to the standard APHA.

Title	Part Number
Dipslide Comparator - App	
DI10 Incubator	56B000701