Lovibond® Colour Measurement Tintometer® Group



FAC

always be reported.



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The valu-

es achieved and the performance tolerances expected will, of course,

FAC

A single glass conformance filter for checking PFXi instruments in-FAC scale. This filter was previously provided as part of the accessories kit for the PFX195/3.

Fiche technique

Glass Standard Valu- Valeurs approximatives : 9

Code: 131043

Filters

Contenu de livraison

FAC 9 filter

Certificate of conformity