

# Lovibond® Water Testing Tintometer® Group



## Chlorine Dioxide LR Strips

For fast determination with chemical reaction



- Hazard free
- Rapid
- Always to hand
- Semi quantitative results

Part Number: 56S000450

### Industry

Industries Others | Municipalities | NGO

### Application

Disaster Management | Disinfection Control

### Chlorine Dioxide LR Strips

Test strips offer quick, clean results without the fuss of reagents or tablets. Their dip and compare method makes them a popular choice for many applications where speed and simplicity are most important.

### Measuring Range

Test Name	Measuring Range
Chlorine Dioxide	0.1 - 1.6 mg/L ClO <sub>2</sub>