

# Lovibond® Colour Measurement Tintometer® Group



## Indice ASTM - EC3000

Code: 135089



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Indice ASTM - EC3000

A single glass filter for checking the EC3000 instrument in ASTM scale.

### Fiche technique

Glass Standard Values	- 4.7
Filters	1

### Contenu de livraison

- Filter
- Certificate of Conformity