

Lovibond® Water Testing

Tintometer® Group



DBNPA Test Kit

For the standard situation in water testing



- For testing cooling water and other industrial applications
- Test equipment with Checkit® Comparator
- Reagents for up to 100 tests
- Easy to use & result evaluation with visual comparison

N. ordine: 56K701190

A fast biocide

DBNPA (dibromonitrile propionamide) is a fast-acting biocide that rapidly degrades to non-hazardous by-products, including ammonia, bromide ions, dibromoacetonitrile and dibromoacetic acid. It therefore also meets strict environmental requirements. It is one of the biocides recommended by many guidelines for the treatment of industrial process water.

Reliably eliminating microorganisms

In water, DBNPA kills microorganisms quickly and acts reliably against bacteria, fungi, yeasts, algae and microbiological slime: it can quickly penetrate the microbial cell membrane, act on certain protein genes and thus lead to cell death.

Test Kit with comparator

The application of the Test Kit is very easy. After pre-filtration of the water sample, a certain number of reagents is added step by step. Further filtration and the addition of reagent tablets follows. The colour change is used to interpret the results with the help of the colour disc in the comparator.

Industria

Altre industrie

Applicazione

Controllo disinfettante

DBNPA Test Kit

The use of DBNPA in cooling water treatment is widespread because it is a fast and reliable biocide. It is used to combat microorganisms that can cause damage to the water system and material. To ensure that it can work efficiently, regular checks of the DBNPA content should be carried out. This test kit is a reliable companion for these routine tasks.

Campo di misura

Test Name

DPNPA

Campo di misura

0 - 6.8 mg/L DBNPA

Fornitura

- Plastic case
- CHECKIT® Comparator with standard optics
- CHECKIT® Disc DBNPA
- 30 tablet reagents DPD No.1 rapid
- 30 tablet reagents DPD No.3 rapid
- 2 Plastic tube with lid, 10 mL
- Instruction manual

Titolo	N. ordine
CHECKIT® Comparator	145000
Cuvetta in plastica, profondità strato 13,5 mm, con coperchio grigio, set da 10	145500
Cuvetta in plastica, profondità strato 13,5 mm, con coperchio grigio, set da 5	145505
Plastic cell, 13,5 mm path length, with grey lid, set of 100	145510
3/150	255020
Siringa, plastica, 20 mL	56A006501

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germania

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
sales@lovibond.uk
www.lovibond.com
Regno Unito

Tintometer China

Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
Cina

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malesia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brasile

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
Stati Uniti d'America

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spagna