

# Lovibond® Colour Measurement Tintometer® Group



## Отдельный фильтр (сертифицированный)

Номер заказа: 109970



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Отдельный фильтр (сертифицированный)

This part number stands for a single filter supplied in other filter sets. Please note that each individual filter has its own part number. Please contact [sales@tintometer.com](mailto:sales@tintometer.com) for the correct number for your filter.

### Технические данные

<b>Glass Standard Values</b>	specify instrument, scale and nominal value when ordering.
<b>Filters</b>	1

### Объем поставки

- 1 filter usually supplied in another set
- Certificate of conformity