

# Lovibond® Colour Measurement Tintometer® Group



## Единицы цветности по шкале ASTM 1

Номер заказа: 134000



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Единицы цветности по шкале ASTM 1

Бутылка 500 мл, сертифицирована согласно ISO/IEC 17025, срок годности 12 месяцев.

Возможность полного отслеживания согласно признанным международным стандартам

Каждая бутылка поставляется с полным комплектом сертификатов, включая паспорт безопасности.