

# Lovibond® Colour Measurement Tintometer® Group



## Цветность по шкале Lovibond® RYBN 2,5R 14,0Y 0,7N (5¼)

Номер заказа: 134120



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Цветность по шкале Lovibond® RYBN 2,5R 14,0Y 0,7N (5¼)

Бутылка 500 мл, сертифицирована для систем контроля качества согласно ISO 9001, срок годности 12 месяцев.

Возможность полного отслеживания согласно признанным международным стандартам

Каждая бутылка поставляется с полным комплектом сертификатов, включая паспорт безопасности.