

# Lovibond® Colour Measurement Tintometer® Group



## Набор фильтров для шкалы Хазена - ЕС 2000

Номер заказа: 135150



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Набор фильтров для шкалы Хазена - ЕС 2000

Glass conformance filters for checking Pt-Co scale on EC 2000 PT-Co.

#### Технические данные

<b>Glass Standard Values</b>	Approximate values - 65 - Zero Baseline
------------------------------	-----------------------------------------------

<b>Filters</b>	2
----------------	---

### Объем поставки

- 65, Zero Baseline filters
- Certificate of Conformity