

Lovibond® Colour Measurement Tintometer® Group



Набор фильтров для оценки меда (шкала Пфунда)

Номер заказа: 139370



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Набор фильтров для оценки меда (шкала Пфунда)

A set of five glass conformance filters for checking PFXi instruments in Pfund honey scale.

Технические данные

Glass Standard Values	Approximate values
	- 15
	- 30
	- 60
	- 85
	- 100

Filters	5
---------	---

Объем поставки

- Pfund 15mm filter
- Pfund 30mm filter
- Pfund 60mm filter
- Pfund 85mm filter
- Pfund 100mm filter
- Certificate of conformity