

Lovibond® Colour Measurement Tintometer® Group



Набор фильтров для шкалы EP желтый

Номер заказа: 139420



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Набор фильтров для шкалы EP желтый

A set of two glass conformance filters for checking PFXi instruments in European Pharmacopoeia Red scale.

Технические данные

Glass Standard - Y2
Values - Y4

Filters 2

Объем поставки

- EP Yellow 2 filter
- EP Yellow 4 filter
- Certificate of Conformity