

Lovibond® Colour Measurement Tintometer® Group



Набор фильтров для шкалы Американской фармакопеи (US Pharmacopoeia)

Номер заказа: 139460



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Набор фильтров для шкалы Американской фармакопеи (US Pharmacopoeia)

A set of three glass conformance filters for checking PFXi instruments in US Pharmacopoeia scale.

Технические данные

Glass Standard Values	G, H, P
Filters	3

Объем поставки

- USP G filter
- USP H filter
- USP P filter
- Certificate of conformity