

Lovibond® Colour Measurement Tintometer® Group



Набор фильтров для шкалы Lovibond® красный, желтый, нейтральный (PFXi 995/950/880)

Номер заказа: 139590



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Набор фильтров для шкалы Lovibond® красный,
желтый, нейтральный (PFXi 995/950/880)

A set of ten glass conformance filters for checking PFXi
instruments (PFXi 995/950/880) in Lovibond RYBN scale.

Технические данные

Glass Standard Values	Approximate Values
	- 0.1R 0.5Y 0.0N
	- 0.5R 1.3Y 0.1N
	- 1.0R 2.6Y 0.1N
	- 1.3R 3.3Y 0.5N
	- 2.0R 5.7Y 0.1N
	- 2.5R 8.6Y 0.1N
	- 3.2R 15.0Y 0.1N
	- 4.0R 22.0Y 0.1N
	- 5.0R 34.0Y 0.1N
	- 7.1R 51.0Y 0.5N

Filters 10

Объем поставки

- Lovibond RYBN 0.1R 0.5Y 0.0N filter
- Lovibond RYBN 0.5R 1.3Y 0.1N filter
- Lovibond RYBN 1.0R 2.6Y 0.1N filter
- Lovibond RYBN 1.3R 3.3Y 0.5N filter
- Lovibond RYBN 2.0R 5.7Y 0.1N filter
- Lovibond RYBN 2.5R 8.6Y 0.1N filter
- Lovibond RYBN 3.2R 15.0Y 0.1N filter
- Lovibond RYBN 4.0R 22.0Y 0.1N filter
- Lovibond RYBN 5.0R 34.0Y 0.1N filter
- Lovibond RYBN 7.1R 51.0Y 0.5N filter
- Certificate of Conformity

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Германия

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
sales@lovibond.uk
www.lovibond.com
ВЕЛИКОБРИТАНИЯ

Tintometer China

Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
Китай

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebu Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Малайзия

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Бразилия

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
США

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
Индия

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Испания

Возможны технические изменения

Отпечатано в Германии

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies