

# Lovibond® Colour Measurement Tintometer® Group



## Набор фильтров для шкалы Lovibond® красный, желтый, синий (PFXi 995/950/880)

Номер заказа: 139620



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Набор фильтров для шкалы Lovibond® красный, желтый, синий (PFXi 995/950/880)

A set of five glass conformance filters for checking PFXi instruments (PFXi 995/950/880) in Lovibond RYBN scale.

#### Технические данные

Glass Standard Values	Approximate Values
	- 1.0R 0.5B
	- 1.3R 0.8B
	- 1.2R 4.3Y 1.2B
	- 2.6R 11.5Y 2.1B
	- 3.4R 15.9Y 2.9B

Filters 5

### Объем поставки

- Lovibond RYBN 1.0R 0.5B filter
- Lovibond RYBN 1.3R 0.8B filter
- Lovibond RYBN 1.2R 4.3Y 1.2B filter
- Lovibond RYBN 2.6R 11.5Y 2.1B filter
- Lovibond RYBN 3.4R 15.9Y 2.9B filter
- Certificate of conformity.