

Lovibond® Colour Measurement Tintometer® Group



Набор фильтров для шкалы Lovibond® красный, желтый, нейтральный (Fx)

Номер заказа: 169610



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Набор фильтров для шкалы Lovibond® красный, желтый, нейтральный (Fx)

A set of five glass conformance filters for checking Model Fx instruments in Lovibond RYBN scale.

Технические данные

Glass Standard Values	Approximate Values
	- 0.1R 0.5Y 0.0N
	- 1.0R 2.6Y 0.1N
	- 2.0R 5.7Y 0.1N
	- 3.2R 15.0Y 0.1N
	- 5.0R 34.0Y 0.1N

Filters 5

Объем поставки

- Lovibond 0.1R 0.5Y 0.0N filter
- Lovibond 1.0R 2.6Y 0.1N filter
- Lovibond 2.0R 5.7Y 0.1N filter
- Lovibond 3.2R 15.0Y 0.1N filter
- Lovibond 5.0R 34.0Y 0.1N filter
- Certificate of conformity