

Lovibond® Colour Measurement Tintometer® Group



ASTM renk endeksi 7

Ürün No: 134030



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

ASTM renk endeksi 7

500 ml şişe, ISO/IEC 17025 uyarınca 12 aylık kullanım ömrüyle sertifikalanmıştır.

Uluslararası kabul görmüş standartlara göre tam takip edilebilirlik

Her şişe, güvenlik bilgi formu dahil tam sertifikalarla birlikte teslim edilir.