

# Lovibond® Colour Measurement Tintometer® Group



## Saybolt renk endeksi +12

Ürün No: 134060



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Saybolt renk endeksi +12

500 ml şişe, ISO/IEC 17025 uyarınca 12 aylık kullanım ömrüyle sertifikalanmıştır.

Uluslararası kabul görmüş standartlara göre tam takip edilebilirlik

Her şişe, güvenlik bilgi formu dahil tam sertifikalarla birlikte teslim edilir.