

# Lovibond® Colour Measurement Tintometer® Group



Lovibond® RYBN renk endeksi 2,5R 14,0Y  
0,7N (5¼)

Ürün No: 134120



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

Lovibond® RYBN renk endeksi 2,5R 14,0Y 0,7N  
(5¼)

500 ml şişe, ISO 9001 uyarınca 12 aylık kullanım ömrüyle QE sistemine göre sertifikalanmıştır.

Uluslararası kabul görmüş standartlara göre tam takip edilebilirlik

Her şişe, güvenlik bilgi formu dahil tam sertifikalarla birlikte teslim edilir.