

Lovibond® Colour Measurement Tintometer® Group



Platin-Kobalt renk endeksi sıvı standardı 15 - 500 ml

Ürün No: 134160



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

Platin-Kobalt renk endeksi sıvı standardı 15 - 500 ml

500 ml şişe, ISO 9001 uyarınca 12 aylık kullanım ömrüyle QE sistemine göre sertifikalanmıştır.

Uluslararası kabul görmüş standartlara göre tam takip edilebilirlik

Her şişe, güvenlik bilgi formu dahil tam sertifikalarla birlikte teslim edilir.