

Lovibond® Colour Measurement Tintometer® Group



Lovibond Red, Yellow, Neutral

Part Number: 139592



Conformance filters are assigned a nominal value and the company always endeavors to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Lovibond Red, Yellow, Neutral

A single glass conformance filter for checking PFXi instruments in Lovibond RYBN scale. This filter is supplied as part of the PFXi 880/F package.

Technical Data

Glass Standard Values	Approximate Value
	- 0.5R
	- 1.3Y
	- 0.1N

Filters	1
---------	---

Delivery Scope

- Lovibond 0.5R, 1.3Y, 0.1N filter
- Certificate of conformity