

Lovibond® Water Testing Tintometer® Group



HARDNESS-total TEST KIT GH-1

For the standard situation in water testing



- Rapid quantitative analysis
- Ideal for boiler, cooling and process water
- Field and laboratory tests possible
- Economical due to inexpensive refill units

Part Number: 2418711

HARDNESS-total TEST KIT GH-1

These test kits are specially developed for testing boiler, cooling and industrial process water. They make use of both colorimetric and titrimetric techniques. Each test kit contains all the necessary chemicals and reagents in liquid or powder form to conduct the tests. The detailed instructions contain a step-by-step explanation of the test procedure. The kits are supplied in a sturdy, compact plastic case. Keenly priced refill reagent packs are available for all test kits.

Measuring Range

Test Name	Measuring Range
Hardness total	1 drop = 1 or 0.5 °dH

Technical Data

Dimensions 90 x 72 x 36 mm

Title	Part Number
Hardness-Total-Reagent	418563