

Lovibond® Water Testing

Tintometer® Group



Legionella Single Syringe EU Test Kit

The World's Fastest *Legionella* Test



Legionella is the causative agent of Legionnaire's disease, a potentially fatal form of pneumonia. Guidelines exist across the world for the control of *Legionella*.

In 2002, the Health and Safety Executive released guidelines for the control of *Legionella*, known as the L8 guidelines. These guidelines detail the key parameters that should be tested on a weekly and monthly basis to minimise the risk of *Legionella* proliferating to harmful levels in open water system.

The Lovibond® *Legionella* "L8" Test Kits contain the recommended tests for both weekly and monthly monitoring.

- Fast, easy and accurate
- Results in as little as 30 minutes!
- Confirms performance of *Legionella* prevention and treatment protocols
- Test Kits available for detecting both waterborne and biofilm *Legionella*

Part Number: 56B006105

Legionella Single Syringe EU Test Kit

The Single Test Kit is for detecting *Legionella* in a water sample. is ideal for users who require the increased test accuracy associated with a filtered sample.

With results in 30-35 minutes, the test allows for immediate understanding and response to the water conditions. Regular use will increase confidence in your *Legionella* control procedures and reduce risk of outbreaks and litigation!

Kit supplies everything required for 1 test.

Limit of Detection
1,000 / 10,000 CFU/l

No. of Tests
1

Delivery Scope

- 1 Test Strip with Dispensing Pipette
- Recovery Buffer Syringe
- Sample Collection Jar
- Single-use filter
- 60 ml syringe
- Visual Score Card
- Instruction Manual