

# Lovibond® Water Testing Tintometer® Group



## Copper Strips 0.05-2 mg/l (25 pc.)

For fast determination with chemical reaction



- Hazard free
- Rapid
- Always to hand
- Semi quantitative results

Part Number: 56S000625

### Industry

Industries Others

### Application

Cooling Water

### Copper Strips 0.05-2 mg/l (25 pc.)

Test strips offer quick, clean results without the fuss of reagents or tablets. Their dip and compare method makes them a popular choice for many applications where speed and simplicity are most important.

### Measuring Range

**Test Name**

Copper-Test Strips

**Measuring Range**

0.05 - 2 mg/L Cu