

# Lovibond® Water Testing

## Tintometer® Group



## Pooltester Chloor LR/pH in doos

The fast assistants for Pool Testing



- Easy to use
- RAPID tablets fast dissolving
- Highest accuracy

Bestelnr.: 151613

### High measuring accuracy

Results are shown with high measurement accuracy. They have a certainty and can counteract if necessary.

### Intuitive

Self-explanatory is the operation of the tester. And if you are ever unsure: A look at the packaging and the intuitively understandable instructions are sufficient to eliminate the ambiguity.

### Ergo mixed

The ergonomically optimized lid ensures a secure grip.

### Industrie

Privézwembad

### Applikation

Behandeling drinkwater | Controle desinfecteermiddel | Controle-zwembadwater | Koelwater | Rampenbeheer | Voedingsmiddelen en dranken

### Pooltester Chloor LR/pH in doos

Een stap naar een duurzame toekomst: de Pooltester CI LR, pH is nu verkrijgbaar in een recyclebare verpakking.

Met de Pooltester kun je gelijktijdig het chloorgehalte en de pH-waarde meten: vul het zwembadwater bij, voeg Rapid-tabletten toe en vergelijk de kleurverandering met de schaalverdeling. Dit maakt het gemakkelijk voor zwembadeigenaren om het wateronderhoud van hun zwembad bij te houden.

### Meetbereik

Test Name	Meetbereik
Chlorine	0.1 - 3.0 mg/L Cl <sub>2</sub>
pH	6.8 - 8.2 pH

# Leveromvang

- Pooltester
- Tablet reagents for 20 tests
- Instruction manual

## **Tintometer GmbH**

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Duitsland

## **The Tintometer Limited**

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
Verenigd Koninkrijk

## **Tintometer China**

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

## **Tintometer South East Asia**

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Maleisië

## **Tintometer Brazil**

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazilië

## **Tintometer Inc.**

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
Verenigde Staten

## **Tintometer India Pvt. Ltd.**

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

## **Tintometer Spain**

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spanje