


# Lovibond® Water Testing

## Tintometer® Group



## Water Safety Kit Basic

Uitgerust voor alle vereisten



Don't wait for solutions in disasters or difficult environments - simply react. With the Lovibond® Water Safety Kit Basic you take everything with you to manage each testing challenge and provide safe potable water. The kit combines microbiological and simple chemical methods of analysis.

Perform microbiological membrane filtration and apply chemical test equipment directly on-site - rapidly & with reliable results.

Its compact design and the optional backpack ensures your mobility in each environment and situation. The Heating and cooling function of the incubator allows you test different sorts of bacteria. Be sure that you don't need any further test kit or test equipment to control the most important parameters you need for safe potable water quality.

- Emergency response kit: Manage each testing challenge
- Portable basic test equipment: Microbiological & chemical analysis in one kit
- Incubator with heating and cooling as standard
- Safe & simple testing according to WHO standards

Bestelnr.: 56K681250

### Industrie

Gemeenten | Ngo | Openbaar zwembad | Overige industrieën | Scheepvaart | Voedingsmiddelen- en drankenindustrie

### Applikation

Behandeling drinkwater | Controle desinfecteermiddel | Koelwater | Rampenbeheer | Voedingsmiddelen en dranken

### Water Safety Kit Basic

Whether disaster or monitoring in difficult environments: You need mobile testing equipment for portable water in each situation. the Water Safety Kit Basic accompanies you with the complete test methods for the most important parameters of both microbiological and physio chemical analysis.

Profit from the advantages of rugged, waterproof cases with compact design and easy to use content.

### User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

## Meetbereik

Test Name	Meetbereik
Aërobe bacteriën (TVC)/colibacteriën	
Aërobe bacteriën (TVC)/E.Coli, coliformen	
Chlorine	0.1 - 3.0 mg/L Cl <sub>2</sub>
Conductivity	0 - 20 mS/cm
pH	6.8 - 8.2 pH
Turbidity	30 - 400 NTU

## Technische gegevens

**Afmetingen** 487 x 229 x 386 mm

**Gewicht met verpakking** 16 kg

Omschrijving	Bestelnr.
Stroomkabelset	136300
Incubator Battery, Lead/Acid	192345
Flesjes	363590
Universal container with spoon and lid, 30 ml	366239
Plastic beker, gegraduateerd, 100 ml, met deksel	385130
Petrischaaltjes, 10 stuks	400855
Troebelheidstest 5-500 JTU	56A010801
Rugzak voor Water Safety Kits	56A01400
Media Pads and Dispenser (Pack of 200)	56A016330
Maatbeker	56A018501
Membraanfiltratie	56A023801
Pincet	56A604201
Turbidity tube (one piece) - 30-400 NTU	56A607101
Membraanfiltratie	56B000710
D004 TTC/MAC Dipslides (Macconkey)	56B010410
D005 TTC/Chromogeen Dipslides	56B010510
KS16-1413US Standaardgeleidbaarheidsoplossing	56L001665
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

## Leveromvang

- DI 20 Incubator
- Incubator battery
- 50 Petri dish for DI 20 incubator
- SD 70 Conductivity Meter
- Marker pen
- 2.5 " Bailer with rope
- Filter flask, 500 ml
- Tweezer
- Pump (35 ml)
- Universal container and spoon, 30 ml
- DB 513 Shaker tube with lid, 100 ml
- 150 EZ Pack filters
- Pipette disposable 1.5 ml (EDT222)
- Turbidity Tube
- KS16 1413 µS Conductivity standard solution (65 ml)
- 200 Media pads and dispenser
- Filtration device
- Lead set IEC C13 (3 pin)
- Membrane Lauryl Sulphate Broth (40 g)
- Chlorine-pH Tester
- 250 DPD No.1 Photometer tablet reagents
- 250 DPD No.3 Evo Photometer tablet reagents
- 250 Phenol Red Photometer tablet reagents
- Instruction manual for Water Safety Kit Basic
- Instruction manual for DI 20 Incubator
- Case

---

**Tintometer GmbH**

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Duitsland

**The Tintometer Limited**

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
Verenigd Koninkrijk

**Tintometer China**

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

**Tintometer South East Asia**

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Maleisië

**Tintometer Brazil**

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazilië

**Tintometer Inc.**

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
Verenigde Staten

**Tintometer India Pvt. Ltd.**

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

**Tintometer Spain**

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spanje