


# Lovibond® Water Testing

## Tintometer® Group



## Closed System Engineers Test Kit

The expert tool for reliable result



These Test Kits for cooling water monitoring are specially designed for engineers and mobile use. It provides the equipment and reagents required to monitor the correct amount and dosing levels of various important parameters. You can carry everything you need with you and test directly onsite. The Kit is provided in a robust carry case and includes all important equipment from measurement instruments to reagents.

- Complete analysis kit including molybdate, iron and conductivity
- Prevents corrosion and scale
- Accurate inhibitor dosage rates
- Reproducible, accurate results

Artigo No: 56K701600

### Indústria

Eletricidade e Energia | Indústria de Alimentos e Bebidas | Indústria-Farmacêutica | Indústria Naval | Indústria Química | Municipalidades | Óleo e Gás | Outras Indústrias

### Aplicação

Água de Caldeira | Água de Refrigeração | Tratamento de Esgotos

### Closed System Engineers Test Kit

Due to its heat absorption capacity, water is the preferred cooling medium in many areas of industry. However, water always contains many dissolved and undissolved substances that can cause trouble even when travelling in closed cooling circuits. Extensive water treatment or chemical post-treatment by experts is always necessary. Problems are mainly caused by deposits, corrosion and microbiological activities. For the experts, a regular check of certain parameters is therefore just as essential as reliable equipment.

## Faixa de Medição

Test Name	Faixa de Medição	Método Químico
Alkalinity M	50 - 2400 mg/L CaCO <sub>3</sub>	
Alumínio T	0 - 0.3 mg/L Al	Eriochrom Cyanine R
Chloride	20 - 12000 mg/L Cl <sup>-</sup>	
Cobre HR T	0 - 5 mg/L Cu	Biquinoline
Dureza, total	5 - 600 mg/L CaCO <sub>3</sub>	
Ferro LR T	0.05 - 1 mg/L Fe	PPST
Glicol	% PEG/MEG	
Molibdato HR T	50 - 500 mg/L MoO <sub>4</sub>	Thioglycolate
Nitrito	10 - 2000 mg/L NaNO <sub>2</sub>	

## Texto do Escopo de Entrega

- Plastic case
- SD40 Multi Pocket Tester
- Checkit Comparator
- 50 - 500 mg/L Molybdate HR+ Checkit Disc
- 0 - 1.0 mg/L Iron LR Checkit Disc
- 0 - 0.3 mg/L Aluminium Checkit Disc
- 0 - 5.0 mg/L Copper HR Free and Total Checkit Disc
- 2 Plastic vials 10 mL with lid
- Aluminium No.1 tablets
- Aluminium No.2 tablets
- Copper No. 1 Test tablets
- Copper No. 2 Test tablets
- Molybdate No.1 tablets
- Molybdate No.2 tablets
- pH Buffer 4 solution
- pH Buffer 7 solution
- pH Buffer 10 solution
- Iron Reagent FE5
- Hardness Total Buffer TH2
- Hardness Total Indicator TH1P
- Hardness LR Titrant TH3
- Hardness HR Titrant TH4
- 1413µS Standard Conductivity Solution
- Alkalinity 4.5 Indicator TA4
- Alkalinity LR Titrant TA3
- Alkalinity HR Titrant PA2/TA2
- Nitrite Indicator N1
- Nitrite Titrant N2
- Chloride Indicator CC1
- Chloride LR Titrant CC2
- Chloride HR Titrant BC2
- 2 Titration jars with caps, plastic, 60 mL
- Syringe plastic, 20 mL
- Wash Bottle, 250 mL
- Syringe Filter Holder 25 mm
- Filter Circle 0.45 µm, 25 mm
- Glycol Refractometer, RHA-218ATC
- Instruction manual

<b>Título</b>	<b>Artigo No</b>
Comparador CHECKIT® Comparador	145000
Seringa, plástico, 20 mL	56A006501
Recipiente de titulação com tampa, plástico, 60 mL	56A006701
Seringa de plástico, 5 mL	56A008501
Frasco de lavagem 250 ml	56A009001
Seringa de plástico 1 ml	56A013501
Círculo de filtragem 0,45 um, 25 mm	56A020050
Tubo de ensaio 5/10 mL + CAP	56A600401
Refratômetro de glicol	56M000401

#### **Tintometer GmbH**

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Alemanha

#### **The Tintometer Limited**

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
Reino Unido

#### **Tintometer China**

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

#### **Tintometer South East Asia**

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malásia

#### **Tintometer Brazil**

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brasil

#### **Tintometer Inc.**

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
EUA

#### **Tintometer India Pvt. Ltd.**

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
Índia

#### **Tintometer Spain**

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Espanha

Technical changes without notice

Printed in Germany

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies