

Lovibond® Water Testing Tintometer® Group



Water Safety Kit Microbiology Duo

Оснащение для любых требований



- Testing coliforms and E. coli bacteria according to WHO standards
- Testing different forms of bacteria simultaneously with 2 incubators
- DI 20 with heating and cooling as standard

Номер заказа: 56K681255

Disaster situation or potable water monitoring in demanding situations: You need certainty and reliability for test results in everywhere. The Water Safety Kit Microbiology Duo contains all the equipment necessary to perform the microbiological analysis component for safe drinking water. Due to the heating and cooling function of the DI 20 incubator you can test different sorts of bacteria - with both incubators simultaneously! Perform filtration with sterilisable equipment including sampler, filter flask & vacuum pump on-site - sterilisation equipment not included. With the optional backpack you can take the kit everywhere with you.

Промышленность

NGO | Другие отрасли

Применение

Подготовка питьевой воды

Water Safety Kit Microbiology Duo

In one hand or on the back: This test kit accompanies you in every situation with the complete microbiological test equipment. Test different sorts of bacteria directly on-site - with both incubators even simultaneously! Be sure of the quality of potable water everywhere, even when there is no laboratory nearby. The DI 20 incubator with its unique capability of heating and cooling as standard ensures reliable results.

User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

Технические данные

Размеры	487 x 229 x 386 mm
Вес с упаковкой	18.75 kg

Объем поставки

- Marker Pen
- 2.5" Bailer with Rope
- Filter Flask, 500ml
- Tweezers
- Pump, 35ml stroke
- 5 Universal Container and spoon, 30 ml
- DB 513 Shaker Tube with lid, 100 ml
- 150 EZ Pack Filters
- 5 Pipette Disposable, 1.5 ml (EDT222)
- 2 DI 20 Incubator
- 200 Media Pads and dispenser
- Incubator Battery Lead/Acid
- Filtration Device
- 50 Petri Dish for DI 20 Incubator
- Lead Set IEC C13, (3 pin)
- Membrane Lauryl Sulphate Broth, 40 g
- Instruction manual for Water Safety Kit Microbiology
- Instruction manual for DI 20 Incubator
- Case

Заголовок	Номер заказа
Набор кабелей питания	136300
Incubator Battery, Lead/Acid	192345
Бутылки	363590
Universal container with spoon and lid, 30 ml	366239
Пластиковый контейнер, градуированный, 100 мл, с крышкой	385130
Чашки Петри, 10 шт.	400855
Рюкзак для Water Safety Kits	56A01400
Media Pads and Dispenser (Pack of 200)	56A016330
Мерные стаканы	56A018501
Мембранная фильтрация	56A023801
Пинцет	56A604201
Мембранная фильтрация	56B000710
D004 пластинки для микробиологического анализа (агар МакКонки) на TTC/МАС	56B010410
D005 пластинки для микробиологического анализа на TTC/пигмент	56B010510
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Германия

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
sales@lovibond.uk
www.lovibond.com
ВЕЛИКОБРИТАНИЯ

Tintometer China

Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
Китай

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Малайзия

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Бразилия

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
США

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
Индия

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Испания

Возможны технические изменения

Отпечатано в Германии

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies