

# Lovibond® Water Testing

## Tintometer® Group



### Polyacrylate Drop Test Kit

Analyse the content of polyacrylate for industrial water conditioning



- Ideal for mobile treatment engineers
- Compact & portable design
- Easy to follow instructions including pictograms
- Complete with all necessary equipment

Номер заказа: 56K700330

The high quality Lovibond® General Drop Test Kits and reagents are available in a wide range of parameter specific options to suit many analytical requirements. Each self-contained Kit is housed in its own rugged carry case and comes complete with all necessary equipment such as instructions, reagents, syringe and test bottle or test cell for 100 tests.

Derived from laboratory burette methods, the Lovibond® drop test range forms an integral part of the analytical program. The reagents are quick to use and give results fast.

#### Промышленность

Другие отрасли

#### Применение

Охлаждающая вода

#### Polyacrylate Drop Test Kit

Lovibond® offers a wide range of single parameter Drop Test Kits for individual test requirements. They are for use in various water testing applications including potable water, process water, industrial boilers & cooling systems, swimming pool and waste water treatment. Drop Tests work by a simple colour change and are derived from laboratory burette methods.

#### Диапазон измерений

Test Name	Диапазон измерений
Polyacrylate	0 - 20 mg/L PAA

## Объем поставки

- Plastic case
- Polyacrylate Buffer A1
- Polyacrylate Precipitant A2
- Polyacrylate Indicator A4
- Polyacrylate HR Titrant A5
- Polyacrylate LR Titrant A6
- Polyacrylate Solvent A3
- Syringe plastic, 20 mL
- Titration jar, glas
- Syringe Filter Holder 25mm
- Filter Papers 25mm 0.45 Micron (50Pk)
- Instruction manual

Заголовок	Номер заказа
Шприц, пластик, 20 мл	56A006501
Стеклянный сосуд для титрования, 50 мл	56A008101
Держатель фильтра 25 мм	56A009101
Фильтровальный диск 0,45 мкм, 25 мм	56A020050

### Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Германия

### The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
ВЕЛИКОБРИТАНИЯ

### Tintometer China

Room 1001, China Life Tower  
16 Chaoyangmenwai Avenue,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 App. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
Китай

### Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Малайзия

### Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Бразилия

### Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
США

### Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
Индия

### Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Испания

Возможны технические изменения

Отпечатано в Германии

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies