

Lovibond® Water Testing Tintometer® Group



Water Treatment Engineers Test Kit (advanced)

Advanced test equipment for water treatment professionals



- All instruments & reagents for daily routine
- Including advanced photometer MD640 for additional parameters
- Early prevention & detection of risks
- Ensure system performance & operating costs

Номер заказа: 56K701500

Water treatment engineers are engaged daily in the control of water quality in closed circuits, evaporative cooling systems and boiler systems. They need reliable equipment and accurate results to rely on. Both have an important impact on decisions on treatment programs and optimal system performance

This Advanced Engineers Kit contains all you need to protect the life of the plant and operating costs. It allows early detection of corrosion, scale and biofouling and the investigation of increased biological activity.

Includes the following parameters: Conductivity, pH, glycol, aluminium, bromine, chlorine, copper, iron, molybdate, phosphate and alkalinity (M, P, OH).

Промышленность

NGO | Другие отрасли | Коммуны | Нефтяная промышленность | Обслуживание плавательных бассейнов | Общественный бассейн | Пищевая промышленность и производство напитков | Поставщики энергии | Судоходство | Фармацевтическая промышленность | Химическая промышленность | Частный бассейн

Применение

Others | Гальванизация | Контроль воды в бассейне | Контроль-дезинфицирующих средств | Котельная вода | Обработка сточных вод | Обработка сырой воды | Охлаждающая вода | Подготовка-питьевой воды

Water Treatment Engineers Test Kit (advanced)

This kit contains pocket meters, checkit comparator, drop tests and MD640 photometer for advanced testing an extended range of control parameters.

Take the complete equipment for direct on-site analysis in practical case in one hand.

Advances analysis of all types of industrial water with visual, electrochemical and photometric methods.

Диапазон измерений

Test Name	Диапазон измерений	Химический метод
Alkalinity (P, M, OH)	50 - 2400 mg/L CaCO ₃	
Chloride	20 - 12000 mg/L Cl ⁻	
Conductivity	0 - 200 mS/cm	
Hardness Calcium	5 - 600 mg/L CaCO ₃	
Hydrogen Peroxide	15 - 500 mg/L H ₂ O ₂	
pH	-2 - 16	
Phosphonate	0 - 20 mg/L HEDP	
PTSA 2P	10 - 400 ppb	Свечение
Sulphite	25 - 150 mg/L Na ₂ SO ₃	
Tannin	50 - 300 mg/L Tannin	
Алюминий Т	0.01 - 0.3 mg/L Al	Эриохромоцианин Р
Бром Т	0.05 - 13 mg/L Br ₂	DPD
Гликоль	% PEG/MEG	
Железо HR L	0.1 - 10 mg/L Fe	Тиогликолят
Железо LR L (A)	0.03 - 2 mg/L Fe	Феррозин / тиогликолят
Жесткость (Да/Нет)	8 - 20 mg/L CaCO ₃	
Жесткость общая	5 - 600 mg/L CaCO ₃	
Медь Т	0.05 - 5 mg/L Cu ^{a)}	Биквинолин
Молибдат Т	1 - 50 mg/L MoO ₄	Тиогликолят
Нитрит	10 - 2000 mg/L NaNO ₂	
Фосфат HR Т	0.33 - 26 mg/L P	Ванадомолибдат
Хлор Т	0.01 - 6.0 mg/L Cl ₂ ^{a)}	DPD

Объем поставки

- Plastic case
- 6 Cuvettes 24 mm round with lids & sealing rings
- SD 335 Multi
- MD640 Photometer
- Brush, 11 cm length
- DPD No.1 tablets BT
- DPD No.3 tablets
- Glycine tablets
- DPD Nitrite tablets bottle
- Aluminium No.1 tablets
- Aluminium No.2 tablets
- Copper No.1 Test tablets
- Copper No.2 Test tablets
- Molybdate No.1 tablets
- Molybdate No.2 tablets
- Phosphate H.R.P1 tablets
- Phosphate H.R.P2 tablets
- Copper/Zinc L.R.Test Tablets
- Hardness Yes/No tablets
- pH Buffer 4 solution
- pH Buffer 7 solution
- pH Buffer 10 solution
- 1413uS Standard Conductivity Solution
- Iron Reagent FE6
- Iron Reagent FE5
- Alkalinity 4.5 Indicator TA4
- Alkalinity LR Titrant TA3
- Alkalinity HR Titrant PA2/TA2
- Alkalinity OH Reagent
- Hardness Total Buffer TH2
- Hardness Total Indicator TH1P
- Hardness LR Titrant TH3
- Hardness HR Titrant TH4
- Chloride LR Titrant CC2
- Chloride HR Titrant BC2
- Tannin Indicator TN1
- Tannin Titrant TN2
- Acetic acid 20%
- Nitrite Titrant N2
- Phosphonate Titrant P5
- Hydrogen Peroxide Buffer HP1
- Phosphonate Indicator P4L
- Nitrite Indicator N1
- Hydrogen Peroxide HR Titrant HP2
- Iron Reagent FE8
- Phosphonate pH Adjuster P3
- Chloride Indicator CC1
- Calcium Hardness Indicator CH1P
- Calcium Hardness Buffer CH2
- Phosphonate Neutraliser P1/2
- Sulphite Indicator S1
- Sulphite Indicator S2
- Syringe plastic, 20 mL

Заголовок	Номер за-каза
LC 12 для SD 320 / SD 325 Соп, до 200 мСм/см	19805040
ValidCheck Хлор 1,5 мг/л	48105510
Шприц, пластик, 20 мл	56A006501
Банка для титрования с крышкой, пластик, 60 мл	56A006701
Пластиковый шприц, 5 мл	56A008501
Пластиковый шприц, 1 мл	56A013501
Пробирка 5/10 мл + CAP	56A600401
Рефрактометр для гликоля	56M000401
SD 335 Multi (базовый блок)	724830

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Германия

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
ВЕЛИКОБРИТАНИЯ

Tintometer China

Room 1001, China Life Tower
16 Chaoyangmenwai Avenue,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 App. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
Китай

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Малайзия

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Бразилия

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
США

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
Индия

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Испания