

# Lovibond® Water Testing

## Tintometer® Group



### Pooltester Phosphate Test Kit, in plastic case, Phosphate 0-1000 ppb

Determination of Phosphate in low concentration levels



- Highly sensitive
- Easy to use
- Futuristic design
- Highest accuracy

Part Number: 157800

#### Highly sensitive

With a 9 cm pathlength, this visual test kit is sensitive enough to detect low levels of phosphate in water.

#### Low cost

This easy-to-use, low cost, testing solution delivers accurate results for measuring the range of 0 - 1000 ppb (= 0 - 1.0 mg/l PO<sub>4</sub>).

#### Industry

Pool Domestic | Pool Services

#### Application

Pool Water Control

#### Pooltester Phosphate Test Kit, in plastic case, Phosphate 0-1000 ppb

Does your pool have an algae problem? Algae growth is often a result of phosphate in the pool water. In order to combat the issue, a phosphate binding material can be added to the water, combined with a routine practice of monitoring phosphate levels in the pool. It is important that the solution you choose can accurately detect the low concentration of

200 µg/l PO<sub>4</sub>, the level at which an algae bloom can occur.

The Lovibond® Phosphate test kit is up to the task!

#### Measuring Range

Test Name	Measuring Range
Phosphate	0 - 1000 ppb

# Delivery Scope

- Color card
- Stirring rod, length 13 cm
- Test vial & stopper
- 30 Phosphate No.1 tablets
- 30 Phosphate No.2 tablets
- Instruction Manual

## Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Germany

## The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
Fax: +44 (0)1980 625412  
sales@lovibond.uk  
www.lovibond.com  
UK

## Tintometer China

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

## Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebu Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malaysia

## Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazil

## Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
USA

## Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

## Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spain