

Lovibond® Water Testing Tintometer® Group



Kit de prueba, Isotiazolinona, 0-7.5 mg/L

For the control of Isothiazoline concentration



- For monitoring of cooling water & other industrial applications
- Works with colour card
- For up to 100 tests
- Carry out the test in a few steps directly on site

Referencia No: 56K701200

Biocide with versatile effect

Isothiazolinone is a biocide for the control of bacteria, fungi and algae, especially in industrial water and cooling water. It works by separating the proteins even at low doses. This inhibits microbial growth and kills the unwanted "inhabitants". Isothiazolinone has good degradability and no side effects. It is safe to use and has good stabilisation properties. Its use is inexpensive.

Easy to use

Professionals and non-professionals can quickly obtain reliable results with this Test Kit. The test works with the colourimetric method: If isothiazolinone is present in the water, the sample turns blue after the addition of certain reagents. The concentration can be read off from the colour chart. The test can be carried out directly on site and provides certainty within minutes - without a laboratory. The easy-to-understand instructions with pictograms ensure that the test can be carried out by any user, regardless of experience.

Many advantages

Many advantages With our Test Kit you can quickly determine how much active isothiazolinone is effective in the water. This helps you avoid incorrect dosing and possible material damage or system blockages.

Industria

Otras industrias

Aplicaciones

Agua de refrigeración | Control de desinfección

Kit de prueba, Isotiazolinona, 0-7.5 mg/L

Isothiazolinone has a broad spectrum of efficacy against microbiological growth in industrial process water, especially in cooling water. It is already effective at very low concentrations. To ensure effectiveness, regular checks of the isothiazolinone concentration are recommended - to avoid costs and possible damage to the material.

The Lovibond® Non-Oxidising Biocide Test Kits also provide a simple solution to ensure bleed concentrations are in line with environmental consents.

Rango de medición

Test Name	Rango de medición
Isothiazoline	0 - 7.5 mg/L

Volumen de suministro

- Plastic case
- Isothiazolinone DK1
- Isothiazolinone DK2
- Isothiazolinone DK3
- Isothiazolinone DK4
- Isothiazolinone DK5
- Stirring rod and powder spoon
- Plastic tube with lid, 10 mL
- Syringe plastic, 20 mL
- 2 Titration jars with cap, plastic, 60 mL
- Instruction manual

Título	Referencia No
Jeringuilla, plástico, 20 mL	56A006501
Varilla agitadora y cucharilla para polvo	56A006601
Recipiente de valoración con tapa, plástico, 60 mL	56A006701
Tubos de plástico con tapa 10 ml	56A006901
Tubo de ensayo 5/10 mL + tapa	56A600401

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Alemania

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
Reino Unido

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malasia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brasil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
Estados Unidos de América

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
España

Nos reservamos el derecho de cambios técnicos

Printed in Germany

Lovibond® and Tintometer® are Trademarks of the Tintometer Group of Companies