

Lovibond® Water Testing Tintometer® Group



AF 129 水平衡

提供 400 多种不同颜色的色盘



- 提供 400 多种不同颜色的色盘
- 有色和混浊样本的补偿
- 保证玻璃滤色片的耐用性
- 集成棱镜

货号: 411290

可用于不同应用

比色仪 2000+ 可以多样化应用：它用于游泳池以及研究中心、大学或一般的饮用水处理领域。

日光灯

为了能够在不利的能见条件下始终安全稳定地测量，我们建议使用（可选）电源或电池供电的日光灯。如此光线条件始终保持不变。无论在哪里测量。

可实行各种分析

对于比色仪 2000+ 提供 400 多种色盘和多种试剂。其组合得到保障。如果您必须对其进行一次更改，完全可以在不打开比色仪的情况下工作。

集成棱镜

集成棱镜在视觉上将色盘的玻璃标准和有色样本带入视野。

行业

NGO | 公共泳池 | 医药行业 | 市政 | 泳池服务

应用

冷却水 | 污水处理 | 泳池水质控制 | 消毒控制 | 锅炉水

AF 129 水平衡

比色仪 2000+ 是一种用于水分析的通用比色系统。操作舒适，精度和重现性无妥协：比色仪 2000+ 可用于诸多应用。游泳池、研究或饮用水处理只是其中几例。

测量范围

Test Name	测量范围	化学方法
值 T	6.8 - 8.4 pH	苯酚红
氯 L	0.2 - 4.0 mg/L Cl ₂	DPD
碱度 m	10 - 500 mg/L CaCO ₃	
钙硬度	10 - 500 mg/L CaCO ₃	

标准配置文本

The Comparator 2000+ Kit is supplied in a robust carrying case containing:

- Comparator 2000+
- pH comparator disc
- Chlorine comparator disc
- 4 x Comparator vials, 10ml, 13mm path length
- Shaker Tube
- Tablet reagents for each parameter
- Crushing / Stirring Rod
- Vial cleaning brush
- 100ml graduated cylinder
- Instruction sheet
- Langelier Index / Balanced Water Calculation Table

标题

货号

带把手的塑料漏斗

471007

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
德国

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
英国

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
中国

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebu Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
马来西亚

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
巴西

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
美国

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
印度

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
西班牙