

Lovibond® Water Testing Tintometer® Group



Water Safety Kit Combined Duo

Extended special test kit for potable water not only for NGOs



- Emergency response kits for disaster situations and water monitoring outback
- Equipment for chemical and microbiological analysis in two cases

货号: 56K681252

You need safe, fast and reliable test results in difficult situations like disasters or potable water monitoring in rural environments without laboratories? You immediate certainty about the quality of potable water to save health and lives? Take all you need with you and combine microbiological methods and simple chemical methods of analysis.

Take the smaller Water Safety Kit Chemical in one hand and the bigger Water Safety Kit Microbiology on your back in the optional backpack. You are able to test different sorts of bacteria on-site due to the heating and cooling function of our DI 20 petri dish incubator.

行业

NGO | 公共泳池 | 其他行业 | 化工行业 | 医药行业 | 家庭泳池 | 市政 | 油气 | 泳池服务 | 海洋工业 | 能源 | 食品和饮料行业

应用

冷却水 | 原水处理 | 污水处理 | 泳池水质控制 | 消毒控制 | 锅炉水 | 饮用水处理

Water Safety Kit Combined Duo

Combine both microbiological and physio chemical analysis of some of the most critical parameters to determine the suitability of water for drinking. Take the extended complete equipment including 2 incubators for bacteria tests with you everywhere you need immediate and reliable test result.

With their rugged, waterproof cases the kits are easy to transport. Even non technical personal is able to use it due to pictorial instructions.

User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

测量范围

Test Name	测量范围	化学方法
Chlorine T	0.01 - 6.0 mg/L Cl ₂ a)	DPD
Conductivity	0 - 20 mS/cm	
pH	0 - 14	
T 硝酸盐	0.08 - 1 mg/L N	锌还原 / NED
T 铵	0.02 - 1 mg/L N	靛酚蓝
好氧菌 (TVC) / 大肠杆菌、大肠菌群		
好氧菌 (TVC) / 大肠菌群		
浊度	0.01 - 1100 NTU	

标准配置文本

- 1 case Water Safety Kit Microbiology Duo
- 1 case Water Safety Kit Chemical

标题	货号
水上安全套装背包	56A01400
Media Pads and Dispenser (Pack of 200)	56A016330
量杯	56A018501
膜式过滤	56A023801
镊子	56A604201
膜式过滤	56B000710
D004 TTC/MAC Dipslide (麦康凯)	56B010410
D005 TTC/呈色法 Dipslide	56B010510
KS3-pH 10 缓冲溶液 65 ml	56L000365
KS6-pH4 缓冲溶液 65 ml	56L000665
KS9-pH 7 缓冲溶液 65 ml	56L000965
KS16-1413US 标准电导率溶液	56L001665
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

标题	货号
电源线组	136300
Incubator Battery, Lead/Acid	192345
浊度标准 T-CAL 套件 (<0.1, 20, 200, 800 NTU)	194150
瓶子	363590
Universal container with spoon and lid, 30 ml	366239
塑料容器, 100 ml 刻度, 带盖	385130
塑料容器, 30 ml 刻度, 带盖	385132
培养皿, 10个。	400855
ValidCheck 氨 1.5 mg/l	48105510
搅拌棒和粉勺	56A006601

Tintometer GmbH
Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
德国

The Tintometer Limited
Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
英国

Tintometer China
9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
中国

Tintometer South East Asia
Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
马来西亚

Tintometer Brazil
Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
巴西

Tintometer Inc.
6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
美国

Tintometer India Pvt. Ltd.
Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
印度

Tintometer Spain
Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
西班牙