

Lovibond® Water Testing Tintometer® Group



Water Safety Kit Microbiology

Mobile microbiological test equipment for demanding situations



- Testing coliforms and E. coli bacteria according to WHO standards
- DI 20 with heating and cooling as standard
- Also available as a duo version

货号: 56K681254

You need certainty and reliability for test results in emergency situations or outback in areas without laboratories. The Water Safety Kit Microbiology contains all the equipment necessary to perform the microbiological analysis component for safe drinking water. Due to the heating and cooling function of the DI 20 incubator you can test different sorts of bacteria. You can perform filtration with sterilisable equipment including sampler, filter flask & vacuum pump on-site - sterilisation equipment not included. A backpack is available as an option for this kit to allow ease of transportation.

行业

NGO | 其他行业

应用

饮用水处理

Water Safety Kit Microbiology

Take all test equipment with you in one hand or on your back in every situation. This kit includes all you need to determine bacteria directly on-site and monitor the quality of potable water even when you don't have any laboratory nearby. The DI 20 incubator with its unique capability of heating and cooling as standard ensures reliable results.

User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

Technical Data

尺寸	487 x 229 x 386 mm
带包装的重量	14.75 kg

标准配置文本

- Marker Pen
- 2.5" Bailer with Rope
- Filter Flask, 500 ml
- Tweezers
- Pump, 35 ml stroke
- 5 Universal Container and spoon, 30 ml
- DB 513 Shaker Tube with lid, 100 ml
- 150 EZ Pack Filters
- 5 Pipette Disposable, 1.5 ml (EDT222)
- DI 20 Incubator
- 200 Media Pads and dispenser
- Incubator Battery Lead/Acid
- Filtration Device
- 50 Petri Dish for DI 20 Incubator
- Lead Set IEC C13, (3 pin)
- Membrane Lauryl Sulphate Broth, 40 g
- Instruction manual for Water Safety Kit Microbiology
- Instruction manual for DI 20 Incubator
- Case

标题	货号
电源线组	136300
Incubator Battery, Lead/Acid	192345
瓶子	363590
Universal container with spoon and lid, 30 ml	366239
塑料容器, 100 ml 刻度, 带盖	385130
培养皿, 10个。	400855
水上安全套装背包	56A01400
Media Pads and Dispenser (Pack of 200)	56A016330
量杯	56A018501
膜式过滤	56A023801
镊子	56A604201
膜式过滤	56B000710
D004 TTC/MAC Dipslide (麦康凯)	56B010410
D005 TTC/呈色法 Dipslide	56B010510
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
德国

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
英国

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
中国

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
马来西亚

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
巴西

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
美国

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
印度

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
西班牙