


# Lovibond® Water Testing

## Tintometer® Group



### Water Safety Kit Combined

Uitgerust voor alle vereisten



You need safe, fast and reliable test results in difficult situations like disasters or potable water monitoring in rural environments without laboratories? You immediate certainty about the quality of potable water to save health and lives? Take all you need with you and combine microbiological methods and simple chemical methods of analysis. Take the smaller Water Safety Kit Chemical in one hand and the bigger Water Safety Kit Microbiology on your back in the optional backpack. You are able to test different sorts of bacteria on-site due to the heating and cooling function of our DI 20 petri dish incubators.

- Emergency response kits for disaster situations and water monitoring outback
- Equipment for chemical and microbiological analysis in two cases
- Also available as a duo version

Bestelnr.: 56K681251

#### Industrie

Chemische industrie | Energieleveranciers | Farmaceutische industrie | Gemeenten | Ngo | Olie-industrie | Openbaar zwembad | Overige industrieën | Privézwembad | Scheepvaart | Voedingsmiddelen- en drankenindustrie | Zwembadservice

#### Applikation

Afvalwaterzuivering | Behandeling drinkwater | Controle desinfectiemiddel | Controle zwembadwater | Ketelwater | Koelwater | Zuivering- vervuild water

#### Water Safety Kit Combined

Combine both microbiological and physio chemical analysis of some of the most critical parameters to determine the suitability of water for drinking. Take the whole equipment with you everywhere you need immediate and reliable test result.

With their rugged, waterproof cases the kits are easy to transport. Even non technical personal is able to use it due to pictorial instructions.

#### User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

## Meetbereik

Test Name	Meetbereik	Chemische methode
Aërobe bacteriën (TVC)/ colibacteriën		
Aërobe bacteriën (TVC)/ E.Coli, coliformen		
Ammonium T	0.02 - 1 mg/L N	Indofenolblauw
Chlorine T	0.01 - 6.0 mg/L Cl <sub>2</sub> <sup>a)</sup>	DPD
Conductivity	0 - 20 mS/cm	
Nitraat T	0.08 - 1 mg/L N	Zinkreductie / NED
pH	0 - 14	
Turbidity	0.01 - 1100 NTU	

## Leveromvang

- 1 case Water Safety Kit Microbiology
- 1 case Water Safety Kit Chemical

Omschrijving	Bestelnr.
Media Pads and Dispenser (Pack of 200)	56A016330
Maatbeker	56A018501
Membraanfiltratie	56A023801
Pincet	56A604201
Membraanfiltratie	56B000710
D004 TTC/MAC Dipslides (Macconkey)	56B010410
D005 TTC/Chromogeen Dipslides	56B010510
KS3-pH 10 bufferoplossing 65 ml	56L000365
KS6-pH 4 bufferoplossing 65 ml	56L000665
KS9-pH 7 bufferoplossing 65 ml	56L000965
KS16-1413US Standaardgeleidbaarheidsoplossing	56L001665
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

Omschrijving	Bestelnr.
Stroomkabelset	136300
Incubator Battery, Lead/Acid	192345
Satz Trübungsstandards T-CAL (<0,1, 20, 200, 800 NTU)	194150
Flesjes	363590
Universal container with spoon and lid, 30 ml	366239
Plastic beker, gradueerd, 100 ml, met deksel	385130
Plastic beker, gradueerd, 30 ml, met deksel	385132
Petrischaaltjes, 10 stuks	400855
ValidCheck Chloor 1,5 mg/l	48105510
Roerstaaf en poederlepel	56A006601
Rugzak voor Water Safety Kits	56A01400

**Tintometer GmbH**  
 Lovibond® Water Testing  
 Schleefstraße 8-12  
 44287 Dortmund  
 Tel.: +49 (0)231/94510-0  
 sales@lovibond.com  
 www.lovibond.com  
 Duitsland

**The Tintometer Limited**  
 Lovibond House  
 Sun Rise Way  
 Amesbury, SP4 7GR  
 Tel.: +44 (0)1980 664800  
 sales@lovibond.uk  
 www.lovibond.com  
 Verenigd Koninkrijk

**Tintometer China**  
 9F, SOHO II C.  
 No.9 Guanghualu,  
 Chaoyang District,  
 Beijing, 100020  
 Customer Care China Tel.: 4009021628  
 Tel.: +86 10 85251111 Ext. 330  
 Fax: +86 10 85251001  
 chinaoffice@tintometer.com  
 www.lovibond.com  
 China

**Tintometer South East Asia**  
 Unit B-3-12, BBT One Boulevard,  
 Lebu Nilam 2, Bandar Bukit Tinggi,  
 Klang, 41200, Selangor D.E  
 Tel.: +60 (0)3 3325 2285/6  
 Fax: +60 (0)3 3325 2287  
 lovibond.asia@tintometer.com  
 www.lovibond.com  
 Maleisië

**Tintometer Brazil**  
 Caixa Postal: 271  
 CEP: 13201-970  
 Jundiaí – SP  
 Tel.: +55 (11) 3230-6410  
 sales@lovibond.us  
 www.lovibond.com.br  
 Brazilië

**Tintometer Inc.**  
 6456 Parkland Drive  
 Sarasota, FL 34243  
 Tel: 941.756.6410  
 Fax: 941.727.9654  
 sales@lovibond.us  
 www.lovibond.us  
 Verenigde Staten

**Tintometer India Pvt. Ltd.**  
 Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
 Sanathnagar Industrial Estate,  
 Hyderabad, 500018  
 Telangana  
 Tel: +91 (0) 40 23883300  
 Toll Free: 1 800 599 3891/ 3892  
 indiaoffice@lovibond.in  
 www.lovibondwater.in  
 India

**Tintometer Spain**  
 Postbox: 24047  
 08080 Barcelona  
 Tel.: +34 661 606 770  
 sales@tintometer.es  
 www.lovibond.com  
 Spanje