

Lovibond® Water Testing

Tintometer® Group



Pooltester Copper-pH-Active oxygen Active-oxygen 0-10 mg/L O₂ pH 6.8-8.2

The fast assistants for Pool Testing



- Easy to use
- Futuristic design
- RAPID tablets fast dissolving
- Highest accuracy

Part Number: 155235

High measuring accuracy

Results are shown with high measurement accuracy. They have a certainty and can counteract if necessary.

Intuitive

Self-explanatory is the operation of the tester. And if you are ever unsure:
A look at the packaging and the intuitively understandable instructions are sufficient to eliminate the ambiguity.

Ergo mixed

The ergonomically optimized lid ensures a secure grip.

Industry

Pool Domestic | Pool Services

Application

Pool Water Control

Pooltester Copper-pH-Active oxygen Active-oxygen 0-10 mg/L O₂ pH 6.8-8.2

The pool tester enables the simultaneous measurement of the water treatment agent and the pH value.

Measuring Range

Test Name	Measuring Range
Active-oxygen	0 - 10 mg/L O ₂
Copper	0.1 - 1 mg/L Cu
pH	6.8 - 8.2 pH

Delivery Scope

- Pooltester in a sturdy plastic box
- Tablet reagents for 20 tests
- Instruction manual

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germany

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
sales@lovibond.uk
www.lovibond.com
UK

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebu Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
USA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain