

# Lovibond® Water Testing

## Tintometer® Group



## Legionella Test Kit, Industrial, 5 tests

The World's Fastest *Legionella* Test



- Conforms to regulatory guidelines
- Weekly and monthly testing
- Chlorine and Bromine options
- Easy-to-use, uncompromised accuracy

Part Number: 56B006101

Legionella is the causative agent of Legionnaire's disease, a potentially fatal form of pneumonia. Guidelines exist across the world for the control of Legionella.

In 2002, the Health and Safety Executive released guidelines for the control of Legionella, known as the L8 guidelines. These guidelines detail the key parameters that should be tested on a weekly and monthly basis to minimize the risk of Legionella proliferating to harmful levels in open water system.

The Lovibond® Legionella "L8" Test Kits contain the recommended tests for both weekly and monthly monitoring.

### Industry

Chemical Industry | Food and Beverage Industry | Industries Others | Marine Industry | Municipalities | NGO | Oil and Gas | Pharmaceutical Industry | Power and Energy

### Legionella Test Kit, Industrial, 5 tests

The *Legionella* System Test Kit is for detecting Legionella in a water sample — it is ideal for on-line and a variety of industrial applications.

This test procedure utilizes hollow fiber filtration modules which attach directly to the provided sample point fittings (both a 1/2" and 3/4" female threaded fittings are provided). The filtered sample increases the sensitivity of the test.

This filtration step requires no pumping — just connect the sample point fitting to the water system (such as a shower head or on-line sample point) and use the system pressure to drive the filtration. The pre-filled syringes contain the required recovery buffer, so there is no need to measure or dilute solutions.

With results in 30-35 minutes (including sample prep-time), the test allows for immediate understanding and response to the water system conditions. Regular use will increase confidence in your *Legionella* control procedures and reduce risk of outbreaks and litigation!

#### Limit of Detection

100 CFU/l

#### No. of Tests

5

# Delivery Scope

- 5 Test Strips with Dispensing Pipettes
- 5 Syringes with Recovery Buffer
- Sample Collection Jar
- 5 single-use filters
- Female threaded sample point fitting
- Instruction Manual
- Carrying Case

## Tintometer GmbH

Lovibond® Water Testing  
Schleefstraße 8-12  
44287 Dortmund  
Tel.: +49 (0)231/94510-0  
sales@lovibond.com  
www.lovibond.com  
Germany

## The Tintometer Limited

Lovibond House  
Sun Rise Way  
Amesbury, SP4 7GR  
Tel.: +44 (0)1980 664800  
sales@lovibond.uk  
www.lovibond.com  
UK

## Tintometer China

9F, SOHO II C.  
No.9 Guanghualu,  
Chaoyang District,  
Beijing, 100020  
Customer Care China Tel.: 4009021628  
Tel.: +86 10 85251111 Ext. 330  
Fax: +86 10 85251001  
chinaoffice@tintometer.com  
www.lovibond.com  
China

## Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,  
Lebuh Nilam 2, Bandar Bukit Tinggi,  
Klang, 41200, Selangor D.E  
Tel.: +60 (0)3 3325 2285/6  
Fax: +60 (0)3 3325 2287  
lovibond.asia@tintometer.com  
www.lovibond.com  
Malaysia

## Tintometer Brazil

Caixa Postal: 271  
CEP: 13201-970  
Jundiaí – SP  
Tel.: +55 (11) 3230-6410  
sales@lovibond.us  
www.lovibond.com.br  
Brazil

## Tintometer Inc.

6456 Parkland Drive  
Sarasota, FL 34243  
Tel: 941.756.6410  
Fax: 941.727.9654  
sales@lovibond.us  
www.lovibond.us  
USA

## Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2<sup>nd</sup>, 3<sup>rd</sup> & 4<sup>th</sup> Floor  
Sanathnagar Industrial Estate,  
Hyderabad, 500018  
Telangana  
Tel: +91 (0) 40 23883300  
Toll Free: 1 800 599 3891/ 3892  
indiaoffice@lovibond.in  
www.lovibondwater.in  
India

## Tintometer Spain

Postbox: 24047  
08080 Barcelona  
Tel.: +34 661 606 770  
sales@tintometer.es  
www.lovibond.com  
Spain