

# Lovibond® Colour Measurement Tintometer® Group



## Lovibond® RYBN Color 2.0R 7.0Y 0.5N (5¼)

Part Number: 134100



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavors to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Lovibond® RYBN Color 2.0R 7.0Y 0.5N (5¼)

500ml Bottle certified under ISO 9001 quality system with a 12-month shelf life.

Full traceability to internationally recognized standards

Each bottle supplied with full certification including MSDS