

# Lovibond® Colour Measurement Tintometer® Group



## Bêta-carotène

Code: 169520



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

### Bêta-carotène

A glass conformance filter to check the Model Fx range of instruments accuracy in measuring beta Carotene.

### Fiche technique

**Glass Standard Values** 45 ppm avec un trajet optique de 10 mm

**Filters** 1

### Contenu de livraison

- 45ppm filter
- Certificate of Conformity